



Urgent Safety Communication

Possibility of stent folding in Perceval Sutureless Heart Valve

Device/ Product Name:	Perceval Sutureless Heart Valve		
Lot numbers/Serials:	Item #: REF, Product Description: - ICV1208: PVS21, Perceval Sutureless Aortic Heart Valve size S - ICV1209: PVS23, Perceval Sutureless Aortic Heart Valve size M - ICV1210: PVS25, Perceval Sutureless Aortic Heart Valve size L - ICV1211: PVS27, Perceval Sutureless Aortic Heart Valve size XL		
Manufacturer:	LivaNova		
Problem:	Possibility of stent folding due to Perceval valve oversizing.		
Recommendation/Actions:	Prevention Decalcification, to avoid uneven surfaces; Correct sizing, using available information in the IFU; and Ballooning, with the recommendation to pour warm sterile saline (at 37°C) while ballooning. Early detection Visual inspection, checking that the Perceval stent is correctly deployed; and Performing an intraoperative echographic evaluation after Perceval implant to confirm correct positioning and verify valve functionality under beating heart conditions. There are no required actions for patients already implanted with Perceval outside of normal monitoring and treatment.		

SG-1810-22-H 11/06/2018





Devices/Products photo:	4 mm	
Authorized Representative	Company name:	cigalah group
Details	Contact Person:	Hammad Alrowaili
	Phone:	+966 11 4603575
		+966 505404345
	Email:	a-naif9@hotmail.com

You **should** be aware of the mentioned risks in the notice and **contact** the Authorized Representative for corrective action.

Healthcare Professionals should **report** any adverse events suspected to be associated with affected devices above (or other Medical Devices) to:

National Center for Medical Devices Reporting.

Medical Devices Sector

Saudi Food and Drug Authority

Postal Address: Saudi Arabia - Saudi Food and Drug Authority (3292)

North Ring Road - Al Nafal Unit (1)

Riyadh 13312 - 6288

Tel: +966 (11) 2038222 Ext: 2904

Fax: +966 (11) 2757245

Or

Saudi Vigilance

https://ade.sfda.gov.sa/Home/Report

For latest published Recalls/Alerts, please visit (NCMDR Website)

Sincerely, NCMDR Team



SG-1810-22-H 11/06/2018